ENGLISH	Customer Information
General Information 3M™ ESPE™ Filtek™ P60 Posterior Restorative, is a visible-light activated, radiopaque, restorative composite. It is designed for use in posterior restorations. The filler in Filtek P60 restorative is microfilled. The inorganic filler loading is 61% by volume (without slush treatment) with a particle size range of 0.01 to 3.5 µm. Filtek P60 restorative contains BIS-GMA, UDMA, and BIS-EMA resins. A dental adhesive is used to permanently bond the restoration to the tooth structure. The restorative is available in a variety of shades. It is packaged in individual syringes.	No person is authorized to provide any information which deviates from the information provided in this instruction sheet. Caution: U.S. Federal Law restricts this device to sale or use on the order of a dental professional.
	Warranty 3M ESPE warrants this product will be free from defects in material and manufacture. 3M ESPE MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for their application. This product is available within the warranty period, you exclusive remedy and 3M ESPE's sole obligation shall be repair or replacement of the 3M ESPE product.
	Limitation of Liability Except where prohibited by law, 3M ESPE will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.
	Deutsch
	Precautionary Information for Patients This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known allergic reactions to acrylic resins. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.
	Precautionary Information for Dental Personnel This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed. 3M ESPE MSDS information can be obtained from www.3MESPE.com or contact your local subsidiary.
	Instructions for Use Preparation 1. Prophy: Teeth should be cleaned with pumice and water to remove surface stains. 2. Shade Selection: Before isolating the tooth, select the appropriate shade(s) of restorative material using the Filtek P60 shade guide. 3. Isolation: A rubber dam is the preferred method of isolation. Cotton rolls plus an evacuator can also be used.

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Directions for Indirect Procedure For Inlays, Onlays Or Veneers 1. Shade Selection: Choose the appropriate shade(s) of Filtek P60 restorative prior to isolation. If the restoration is of sufficient depth, use of an opaque shade is recommended. 2. Preparation: Prepare the tooth. 3. Impressioning: After preparation is complete, make an impression of the prepared tooth by following the manufacturer's instructions of the impressioning material chosen. An impressioning material, such as manufactured by 3M ESPE, may be used. 4. Laboratory Procedure 2.1 Pour the impression of the preparation with the stone. Place pins at the preparation site at this time if a "triple tray" type of impression was used. 2.2 Separate the cast from the impression after 45 to 60 minutes. Place pins in die and base the cast as for a typical crown and bridge procedure. Mount or articulate the cast to its counter model to an adequate articulator. 2.3 If a second impression was not sent, pour a second cast using the same impression registration. This is to be used as a working model. 2.4 Section out the preparation with a laboratory saw and trim away excess or, expose the margins so they can be easily worked. Mark the margins with a red pencil if needed. Add a spacer at this time if one is being used. 2.5 Soak the die in water, then with a brush, apply a very thin coat of separating medium to the preparation, let it dry somewhat, then add another thin layer. 2.6 Add the first course of composite to the floor of the preparation, stay short of the margins, light cure for 20 seconds. 2.7 Add the second third of composite. Allow for the last third (incisal) to include the contact areas, light cure for 20 seconds. 2.8 Place the die back into the articulated arch, add the last third of composite to the occlusal surface. Overfill very slightly mesially, distally, and occlusally. This will allow for the mesiodistal contacts and the proper occlusal contact when the opposing arch is brought into occlusion with the uncured initial composite. Light cure for only ten seconds, then remove the die to prevent adhering to adjacent surfaces. Finish the curing process. 2.9 With the occlusal contacts already established, begin removing the excess composite from around the points of contact. Develop the inclines and ridges as per remaining occlusal anatomy. 2.10 Care must be taken when removing the prosthesis from the die. Break off small amounts of the die from around the restoration, the die stop should breakaway freely from the cured restoration, until all of the restoration is recovered. 2.11 Using the master die, check the restoration for flash, undercuts, and fit. Adjust as necessary, then polish. 3. Dental Operatory Procedure 3.1 Roughen the interior surfaces of the indirect restoration. 3.2 Cementation: Cement the prothesis using a 3M ESPE resin cement system by following manufacturer's instructions.	Directions for Indirect Procedure For Inlays, Onlays Or Veneers 1. Shade selection: Choose the appropriate shade(s) of Filtek P60 restorative prior to isolation. If the restoration is of sufficient depth, use of an opaque shade is recommended. 2. Preparation: Prepare the tooth. 3. Impressioning: After preparation is complete, make an impression of the prepared tooth by following the manufacturer's instructions of the impressioning material chosen. An impressioning material, such as manufactured by 3M ESPE, may be used. 4. Laboratory Procedure 2.1 Pour the impression of the preparation with the stone. Place pins at the preparation site at this time if a "triple tray" type of impression was used. 2.2 Separate the cast from the impression after 45 to 60 minutes. Place pins in die and base the cast as for a typical crown and bridge procedure. Mount or articulate the cast to its counter model to an adequate articulator. 2.3 If a second impression was not sent, pour a second cast using the same impression registration. This is to be used as a working model. 2.4 Section out the preparation with a laboratory saw and trim away excess or, expose the margins so they can be easily worked. Mark the margins with a red pencil if needed. Add a spacer at this time if one is being used. 2.5 Soak the die in water, then with a brush, apply a very thin coat of separating medium to the preparation, let it dry somewhat, then add another thin layer. 2.6 Add the first course of composite to the floor of the preparation, stay short of the margins, light cure for 20 seconds. 2.7 Add the second third of composite. Allow for the last third (incisal) to include the contact areas, light cure for 20 seconds. 2.8 Place the die back into the articulated arch, add the last third of composite to the occlusal surface. Overfill very slightly mesially, distally, and occlusally. This will allow for the mesiodistal contacts and the proper occlusal contact when the opposing arch is brought into occlusion with the uncured initial composite. Light cure for only ten seconds, then remove the die to prevent adhering to adjacent surfaces. Finish the curing process. 2.9 With the occlusal contacts already established, begin removing the excess composite from around the points of contact. Develop the inclines and ridges as per remaining occlusal anatomy. 2.10 Care must be taken when removing the prosthesis from the die. Break off small amounts of the die from around the restoration, the die stop should breakaway freely from the cured restoration, until all of the restoration is recovered. 2.11 Using the master die, check the restoration for flash, undercuts, and fit. Adjust as necessary, then polish. 3. Dental Operatory Procedure 3.1 Roughen the interior surfaces of the indirect restoration. 3.2 Cementation: Cement the prothesis using a 3M ESPE resin cement system by following manufacturer's instructions.
---	---

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